10 UNITS INSTRUMENT PACK QREM® CYTOKINE

REF. M03-001

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10 UNITS INSTRUMENT PACK QREM® CYTOKINE

REF. M03-001

The product with reference M03-001 is a procedure pack of class IIa and IIb medical devices intended to be used together with the KIT QREM® CYTOKINE, to obtain platelet byproducts such as growth factors, cytokines and biological mediators from a small sample of autologous peripheral blood.

Content

The procedure pack consists of several sterile products: syringes, needles, caps and a box with ampoules of 3.8% sodium citrate (3.13% Ph.Eur respectively 3.8% Erg-B6) that allow up to 10 platelet byproducts to be obtained by being used together with the KIT QREM® CYTOKINE.

#	PRODUCT	UN
1	BD Plastipak™ Luer-Lok™ Syringe 20ml REF.300629 (<i>Jeringa/ Seringa/ Siringa/ Seringue/ Spritzen/ Syringe</i>)	30
2	BD Microlance™ 3 Hypodermic Needle 21G x 1" (0.8 x 25 mm) REF.301156 BD (<i>Aguja/ Agulha/ Ago/ Aiguilles/ Nadel/ Needle</i>)	20
3	PPS Natrium Citricum 3,13 % REF.50200011 (<i>Citrato de Sodio/ Citrato de Sódio/ Citrato di Sodio/ Citrate de Sodium/ Natriumcitrat/ Sodium citrat</i>)	10
4	Venofix Safety 21G 19mm REF.4056506-01 (<i>Aguja/ Agulha/ Ago/ Aiguilles/ Nadel/ Needle</i>)	10
5	Flexible Syringe Cap REF.418200 (<i>Tapón/ Tampa/ Cappuccio/ Capuchon/ Deckel/ Cap</i>)	30

This procedure pack is intended for use only by healthcare professionals in a medical center. The procedure pack contains several sterile devices and must be properly stored and handled under aseptic conditions to avoid any infection or loss of sterility. Once the components have been used, they can cause illness and must be disposed of according to national regulations and those of each center. Reuse may cause infection and other illnesses and/or injuries. Do not autoclave. Do not use any component of the

pack if its package is open or damaged. Keep the procedure pack in a dry place.

Instructions for Use

IMPORTANT: Read these instructions together with the instructions for the KIT QREM® CYTOKINE and the QREM® DEVICE.

Open an ampoule of 3.8% sodium citrate and, using a needle, load 2ml into one of the 20ml syringes. Using the same syringe, carefully draw 18ml of patient's own blood, gently mix for proper homogenization, cover with a cap to avoid contamination and set aside.

Take another syringe and carefully draw 20ml of patient's own blood, immediately connect the syringe to valve A of the KIT QREM® CYTOKINE and gently empty the blood into the container following the instructions of the KIT QREM® CYTOKINE. Remove the rigid adapter from A valve and discard.

Next, connect the syringe previously loaded with the anticoagulant (sodium citrate) to valve F of the KIT QREM® CYTOKINE and gently empty the blood into the corresponding container. Remove the extension from F and discard.

Place the KIT QREM® CYTOKINE inside the QREM® DEVICE and put it into operation following the QREM® DEVICE instructions.

Once the process is finished, extract the KIT QREM® CYTOKINE and screw another syringe to the valve A. Load the syringe by lowering its body until it is full. Keep in mind that the volume obtained ranges between 8 and 12 ml and that the red blood cells located in the lower part of the container must not be aspirated. Once the syringe is filled, connect a needle to perform the infiltration.

SYSTEM/PROCEDURE PACK PRODUCER

Tecnologia Regenerativa Qrem S.L.

Carrer Sant Pere d'Abanto, 4, 3º 1ª

08014 Barcelona (Spain)

Licence nº 7096-PS

IFU-M03-001-EN (ver: September 6, 2022)